



Portable Massage Table

User Manual

[Revision 3.0 November 2016]

RETAIN THIS MANUAL FOR FUTURE REFERENCE
PLEASE READ THIS MANUAL CAREFULLY BEFORE USE

Table of Contents

Safety..... 3

Parts Identification..... 3


Setting Up the Table4

Closing and Transporting the Table.....5

Safety

- If the frame/legs/table structure is damaged in anyway, do not use the table. Contact an authorized service center or technician.
- Maintain the table periodically. Make sure that all screws are secure and the straps are not deteriorated or frayed. Make sure that there are no cracks in the frame.
- **Maximum weight limit = 250kg.** Do not load the table beyond the maximum weight limit. This could result in serious injury to the user.
- Do not use alcohol-based cleaners when cleaning the vinyl, as this may damage the material. Use warm, soapy water and wipe dry.
- Make sure that when using the table, all four legs are set to the same height. Setting the legs at different heights may make the table unstable.
- Do not allow children to assemble or disassemble the massage table unless under adult supervision.




Parts Identification





No.	Name	No.	Name
1	Table (shown extended)	5	Headrest
2	Armrest (2)	6	Headrest Cushion
3	Arm Sling	7	Table Cover (2)
4	Face Hole Pad		

Setting Up the Table

1. Remove the table from the carry bag and place it on the floor so it is sitting on the rubber feet.
2. Release the 2 outer catches (A) locking the table in the folded position.
3. Open the table (B) – the legs will extend as the table opens. Ensure that the table is fully open and the legs are fully extended.
4. Remove all accessories (armrests, face cushion etc) from the underside of the table.
5. Turn the table over (C) so it is standing securely on the legs.
6. Adjust the table height, if required. To do this, push the leg locating pins (D) in so that the lower leg tube can slide up/down in the upper tube. Slide the leg – at each set of location holes, the leg will lock into position. Adjust as required and ensure that all legs are adjusted to the same height.
7. Attach the armrests (E), if required. To do this, slide the armrest rails



8. Attach the headrest (F), if required. To do this, slide the headrest rails into the holes in the end of the table. To adjust the angle and height of the headrest, rotate the locking lever (G) to the "unlock" position, adjust the headrest as required, then rotate the lever to the "lock" position. Attach the headrest cushion to the headrest.
9. Attach the arm sling (H), if required. To do this, loop the strap over the headrest tube, then thread the strap through the buckle as shown. To adjust, loosen the strap in the buckle and pull through as required until the required length is reached, then tighten the strap through the buckle.
10. To raise the backrest (I) (if required), release the middle catch on the end of the table and lift the backrest to the required position. Ensure that the backrest is secure. To lower the backrest, lift it all the way up to release the locking mechanism, then lower it and secure it with the catch.



Closing and Transporting the Table

1. Disassemble all accessories from the table.
2. Lay the table on its side, with the rubber feet on the ground.
3. Place all accessories back in position on the underside of the table and secure using the elastic straps.
4. Gently fold the table up and lock all catches.
5. Place the table into the carry bag and zip it closed.



1) Some experts believe that the incorrect or prolonged use of almost any product may cause serious injury or death. To help reduce your risk of serious injury or death, refer to the information below. For more information, see www.datastreamserver.com/safety

Consult all documentation, packaging and product labelling before use. Note that some products feature documentation available online. It is recommended to print and retain the documentation.

Before each use, check the product for loose/broken/damaged/missing parts, wear or leaks (if applicable). Never use a product with loose/broken/damaged/missing parts, wear or leaks.

Products must be inspected and serviced (if applicable) by a qualified technician every 6 months. This is based on average residential use by persons of average size and strength, and on a property of average metropolitan size. Use beyond these recommendations may require more frequent inspections/servicing.

Ensure that all users of the product have completed a suitable industry recognised training course before being allowed access to the product.

The product has been supplied by a general merchandise retailer that may not be familiar with your specific application or description of application. Be sure to attain third-party approval from a qualified specialist for your application before use, regardless of any assurances from the retailer or its representatives.

This product is not intended for use where fail-safe operation is required. As with any product (for example, automobile, computer, toaster), there is the possibility of technical issues that may require the repair or replacement of parts, or the product itself. If the possibility of such failure and the associated time it may take to rectify could in any way inconvenience the user, business or employee, or financially affect the user, business or employee, then the product is not suitable for your requirements. This product is not intended for use where incorrect operation or a failure of any kind, including but not limited to, a condition requiring product return, replacement, parts replacement or service by a technician may cause financial loss, loss of employee time or an inconvenience requiring compensation.

If this product has been purchased in error when considering the information presented here, contact the retailer directly for details of their returns policy, if required.

